

FISH & SEAFOOD HISTAMINE





SYMMETRIC HISTAMINE

S9030

Lateral Flow Test Kit for the quantitative determination of Histamine in Fresh or Thawed Frozen Raw Fish Samples and Canned fish.

- ✓ AOAC PTM Certification
- ✓ Water extraction
- ✓ Fast & Simple - 3 mins
- ✓ No acylation process required

QUANTITATIVE RAPID TEST FOR HISTAMINE



BIO-SHIELD HISTAMINE

B6048/B6096/B60576

ELISA Test Kit for the quantitative determination of Histamine in Fresh or Thawed Frozen Raw Fish Samples, Canned fish and fish meal.

- ✓ AOAC PTM Certification
- ✓ No acylation process required

FAST ELISA PROCEDURE 20 MINS



1-STANDARD HISTAMINE

D4648

The test provides a cost-effective, and reliable solution to the analyst, including **only the standard 1** (zero standard) in the kit.

