



# STANDARD MYCOTOXINS

---

## AN OPTIMAL ELISA SYSTEM FOR MYCOTOXIN ANALYSIS

---



Requires only one standard to produce accurate results

---

One common extraction for Aflatoxin B1, Total, DON, ZON, Fumonisin



---

**IN JUST 5 MINS**

The most cost-effective  
ELISA method

# PRECISE AND STABLE LIKE A REGULAR ELISA SYSTEM

---



---

1-Standard ELISA series is an innovative ELISA method using a precalibrated curve for the quantification of mycotoxins.

---

The precalibrated standard curve of the 1-Standard method was run multiple times in ProGnosis Biotech facilities and following a thorough examination and validation it is unique for every lot. The user can download the precalibrated curve on Prognosis Data reader.

---