



## 1-STANDARD AFLATOXIN B1

Cat.Number: D4948

### Ring Test Results

Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score
A.I.A. Ring Test Lotto RTB1 240320 Aflatoxin B1 in Corn March 2020	<limit of detection	<limit of detection	1.267
	22.00	35.54	
	13.50	15.85	
	6.80	7.82	
A.I.A. Ring Test Lotto RTB1 290920 Aflatoxin B1 in Corn September 2020	<limit of detection	<0.11	-0.063
	6.50	6.47	
	13.20	13.35	
	22.10	21.71	