



3PR^{mini}

A lightweight, portable reader for all lateral flow formats, cassettes and strips.

Developed by **ProGnosis Biotech**, it delivers quantitative results in a few minutes.



IDEAL FOR FIELD ANALYSIS



INTEGRATION WITH
S-FLOW SOFTWARE



BLUETOOTH & USB
CONNECTIVITY



SHARABLE RESULTS
VIA ANDROID
MOBILE APP



STAND ALONE
BATTERY



Expand Your Testing Capacity

3PR^{quad}

3PR^{duo}

Standalone reader systems that integrate multiple 3PR Mini readers into one platform. Each reader operates independently, while all connected devices and results are managed from a single screen.

Developed by **ProGnosis Biotech**, they combine efficient walk-away testing with a minimal, modern design.



3PR DUO

Up to **6 test strips** and **2 cassettes** simultaneously



3PR QUAD

Up to **12 test strips** and **4 cassettes** simultaneously



**Multiple
walk-away results**



**One single
screen**



**Independent
readings**

Specifications

- ✓ **Standalone system**, no PC required
- ✓ **One** USB-C cable
- ✓ Integration with **S-Flow MultiView** software