


# KONTROLA KVALITETA HRANE

# AFLATOKSIN M1


# SYMMETRIC M1





# SYMMETRIC AFLATOXIN M1

## KVANTITATIVNI BRZI TEST ZA DETEKCIJU AFLATOKSINA M1 U MLEKU



- Revoluciona Symmetric tehnologija
- Za kravlje, kozije, ovčije i UHT mleko
- Precizan, CV <8%
- Opseg kvantifikacije: 8-150 ppt Aflatoksina M1
- Brz & automatizovan proces
- Softver najnovije generacije - jednostavno deljenje rezultata

Validation of LF Symmetric M1  
& Elisa Bio-Shield M1 ES

**ILVO**

Flanders Research Institute for  
Agriculture, Fisheries and Food