

# MILK & DAIRY ANTIBIOTICS



# A BINARY ANSWER IN THE SHORTEST POSSIBLE TIME



Test cold milk right from the truck, refrigerator or bulk tank

Option for visual interpretation

No need for extra disposables



## QUANTUM BT-CC

## QUANTUM BT-CC

W1530/W1560

Detection of  $\beta$ -lactams, Tetracyclines, Ceftiofur and Cefalexin in cow, sheep and goat milk

■ Final Results in 5 min

## QUANTUM BL

W1430/W1460

Detection of  $\beta$ -lactams in fresh cow milk

■ Final Results in 2 min

W1030/W1060

Detection of  $\beta$ -lactams, including Cefalexin, and Tetracyclines in cow, sheep and goat milk

■ Final Results in 5 min

■ Validated by: ILVO | LILC | PIWET



## QUANTUM BT-S

W1130/W1160

Detection of  $\beta$ -lactams, Tetracyclines and Sulfonamides in cow, sheep and goat milk

■ Final Results in 5 min

■ Validated by: ILVO




## QUANTUM BT-SQ

W1230/W1260

Detection of  $\beta$ -lactams, Tetracyclines, Sulfonamides and Quinolones in cow, sheep and goat milk

■ Final Results in 5 min

■ Validated by: ILVO



## QUANTUM BT-SC

W1330/W1360

Detection of  $\beta$ -lactams, Tetracyclines, Streptomycin and Chloramphenicol in cow, sheep and goat milk


■ Final Results in 5 min



## FAST AND ACCURATE RAPID TESTS

Option for visual interpretation

Incubation at 40°C



## SYMMETRIC B+T

S1048/S1096

Detection of  $\beta$ -lactams and Tetracyclines in cow, goat and sheep milk

■ Dipstick Format

■ Validated by: ILVO | LILC | PIWET



## SYMMETRIC BTS

S1248/S1096

Detection of  $\beta$ -lactams, Tetracyclines and Sulfonamides in cow, goat and sheep milk

■ Dipstick Format

■ Validated by: ILVO





**FROM FARM TO CLOUD**


**Get results  
on the go**  
with our  
**Android  
Mobile App**



3 drops of cold milk



No incubation needed



Shareable Results & GPS location

