

# MILK & DAIRY ANTIBIOTICS



# A BINARY ANSWER IN THE SHORTEST POSSIBLE TIME



Test cold milk right from the truck, refrigerator or bulk tank

Option for visual interpretation

No need for extra disposables



Easy-to-use Cassette Format



## QUANTUM BT-CEF

W1030/W1060/W10120

Detection of  $\beta$ -lactams, including **Cefalexin** and **Tetracyclines** in raw cow, sheep and goat milk, pasteurized milk, milk powder and whey

- Final Results in 5 min
- Validated by: ILVO | LILC | PIWET



## QUANTUM BT-SQ

W1230/W1260/W12120

Detection of  $\beta$ -lactams, **Tetracyclines**, **Sulfonamides** and **Quinolones** in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Final Results in 5 min
- Validated by: ILVO



## QUANTUM BT-S

W1130/W1160/W11120

Detection of  $\beta$ -Lactams, **Tetracyclines** and **Sulfonamides** in raw cow, sheep and goat milk, pasteurized milk, milk powder and whey

- Final Results in 5 min
- Validated by: ILVO | PIWET



## QUANTUM BT-SC

W1330/W1360/W13120

Detection of  $\beta$ -Lactams, **Tetracyclines**, **Streptomycin** and **Chloramphenicol** in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Final Results in 5 min



## QUANTUM BT-CC

W1530/W1560/W15120

Detection of  $\beta$ -lactams, **Tetracyclines**, **Ceftiofur** and **Cefalexin** in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Final Results in 5 min
- Validated by: ILVO | PIWET



Option for visual interpretation

Incubation at 40°C

## SYMMETRIC B+T

S1048/S1096

Detection of  $\beta$ -lactams and **Tetracyclines** in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Dipstick Format
- Validated by: ILVO | LILC | PIWET



## SYMMETRIC BTSC

S1448/S1496

Detection of  $\beta$ -lactams, **Tetracyclines**, **Streptomycin** and **Chloramphenicol**, in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Dipstick Format
- Compliant with Russian MRLs (GOST 34107-2017)




## QUANTUM BL

W1430/W1460/W14120

Detection of  $\beta$ -lactams antibiotics in raw cow and goat milk

- Final Results in 2 min
- Validated by: ILVO



## SYMMETRIC BTS

S1248/S1096

Detection of  $\beta$ -lactams, **Tetracyclines** and **Sulfonamides** in raw cow, sheep and goat milk, pasteurized milk, milk powder and whey

- Dipstick Format
- Validated by: ILVO




## SYMMETRIC FLORFENICOL

S1548/S1596

Detection of **Florfenicol** in raw cow, sheep, goat milk and whey

- Dipstick Format





FROM FARM TO CLOUD

# Android Mobile App

Know Instantly.  
Share Immediately.



Only few drops of cold milk

No incubation needed

Shareable Results & GPS location