

# MILK & DAIRY ANTIBIOTICS



# A BINARY ANSWER IN THE SHORTEST POSSIBLE TIME



Test cold milk right from the truck, refrigerator or bulk tank

Option for visual interpretation

No need for extra disposables



Easy-to-use Cassette Format



## QUANTUM BT-CEF

W1030/W1060/W10120

Detection of  $\beta$ -lactams, including Cefalexin and Tetracyclines in raw cow, sheep and goat milk, pasteurized milk, milk powder and whey

- Final Results in 5 min
- Validated by: ILVO | LILC | PIWET

ILVO Negative Results in: 90 sec

## QUANTUM BT-S

W1130/W1160/W11120

Detection of  $\beta$ -Lactams, Tetracyclines and Sulfonamides in raw cow, sheep and goat milk, pasteurized milk, milk powder and whey

- Final Results in 5 min
- Validated by: ILVO | PIWET

ILVO Negative Results in: 90 sec

## QUANTUM BT-SQ

W1230/W1260/W12120

Detection of  $\beta$ -lactams, Tetracyclines, Sulfonamides and Quinolones in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Final Results in 5 min
- Validated by: ILVO

ILVO Negative Results in: 02 min

## QUANTUM BT-SC

W1330/W1360/W13120

Detection of  $\beta$ -Lactams, Tetracyclines, Streptomycin and Chloramphenicol in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Final Results in 5 min

Negative Results in: 02 min

## QUANTUM BT-CC

W1530/W1560/W15120

Detection of  $\beta$ -lactams, Tetracyclines, Ceftiofur and Cefalexin in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Final Results in 5 min
- Validated by: ILVO | PIWET

ILVO Negative Results in: 02 min

## QUANTUM BL

W1430/W1460/W14120

Detection of  $\beta$ -lactams antibiotics in raw cow and goat milk


- Final Results in 2 min
- Validated by: ILVO

ILVO Negative Results in: 60 sec



Option for visual interpretation

Incubation at 40°C



## SYMMETRIC B+T

S1048/S1096

Detection of  $\beta$ -lactams and Tetracyclines in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Dipstick Format
- Validated by: ILVO | LILC | PIWET

ILVO 5 min

## SYMMETRIC BTS

S1248/S1296 S 1296

Detection of  $\beta$ -lactams, Tetracyclines and Sulfonamides in raw cow, sheep and goat milk, pasteurized milk, milk powder and whey

- Dipstick Format
- Validated by: ILVO

ILVO 5 min

## SYMMETRIC BTSC

S1448/S1496

Detection of  $\beta$ -lactams, Tetracyclines, Streptomycin and Chloramphenicol, in raw cow, sheep and goat milk, pasteurized milk and milk powder

- Dipstick Format
- Compliant with Russian MRLs (TP TC 033/2013)

10 min

## SYMMETRIC FLORFENICOL

S1548/S1596

Detection of Florfenicol in raw cow, sheep, goat milk and whey

- Dipstick Format


5 min



FROM FARM TO CLOUD

# Android Mobile App

Know Instantly.  
Share Immediately.



Only few drops of cold milk

No incubation needed

Shareable Results & GPS location