



**SYMMETRIC  
B+T**

# RAPID TEST FOR THE DETECTION OF **β-LACTAMS & TETRACYCLINES**

The latest addition at the **SYMMETRIC SERIES**  
for the detection of antibiotic residues in milk

## ADVANTAGES



**1-Step** procedure



Rapid results in **5 min**



Visual interpretation




Incubator compatibility (40°C)

## 5-MIN PROCEDURE




**1** mix



**2** dip



**3** incubate



**4** read



**PROGNOSIS  
BIOTECH**

ACCURACY THROUGH INNOVATION

[www.prognosis-biotech.com](http://www.prognosis-biotech.com)