

RAPID TEST

FLU_COVID


V1601/V1605/V1620

BE SURE BE SAFE


Rapid Test FLU_COVID is a qualitative, lateral flow test designed to detect the presence or absence SARS-CoV-2 and Influenza virus antigens in nasal or nasopharyngeal swab specimens. It differentiates between Covid-19 (SARS-CoV-2) and the Flu (Influenza, Types A & B).

- One nasal swab for 2 diagnoses (FLU A/FLU B/COVID-19)
- Same protocol of use with the classic rapid test for COVID-19
- Chemical-free sterilised swab
- 4 lines cassette

TEST PROCEDURE



RESULT INTERPRETATION



**RAPID TEST
FLU_COVID**



www.prognosis-biotech.com

✉ info@prognosis-biotech.com ☎ +30 2410 623922

📍 Farsalon 153 | 41335 Larissa, Greece

**PROGNOSIS
BIOTECH**