

KONTROLA KVALITETA HRANE

AFLATOKSIN M1

QUANTUM M1



KVANTITATIVAN BRZI TEST

Za detekciju Aflatoksina M1
u kravljem, kozijem i ovčijem mleku

- **Nije potrebna inkubacija**
- Jednostavna procedura, samo 3 kapi
- Quantum očitavanje 8 min, konačan rezultat 15 min
- Veoma senzitivna, opseg 0-150 ppt

**SA FARME
NA KLAUD**


**BRZI REZULTATI
UZ POMOĆ ANDROID
APLIKACIJE**



3 kapi hladnog mleka



Inkubacija nije potrebna



Kvantitativni rezultati & GPS lokacija

