



**PROGNOSIS
BIOTECH**

RAPID TEST FLU_COVID


V1601/V1605/V1620

BE SURE BE SAFE


Rapid Test FLU_COVID is a qualitative, lateral flow test designed to detect SARS-CoV-2 and/or Influenza A / B virus antigens in nasal or nasopharyngeal swab specimens directly collected. It differentiates between Covid-19 (SARS-CoV-2) and the Flu (Influenza, Types A & B).

- One nasal swab -> 2 diseases (FLU A/FLU B/ COVID-19)
- Same protocol of use with the classic rapid test for COVID-19
- Chemical-free sterile swab
- 4 lines cassette

TEST PROCEDURE



RESULT INTERPRETATION



RAPID TEST
FLU_COVID



www.prognosis-biotech.com

✉ exports@prognosis-biotech.com ☎ +30 2410 623922

📍 Farsalon 153 | 41335 Larissa, Greece

PROGNOSIS
BIOTECH