



PROGNOSIS
BIOTECH

RAPID TEST

STREP A


V2001/V2005/V2020

YOUR RELIABLE CHOICE FOR GROUP A STREPTOCOCCAL DETECTION


- Visual interpretation of the results in **5 mins**
- Simple method, no need for special equipment
- EU regulatory requirements compliant

Rapid Test STREP A is a qualitative lateral flow test designed for the detection of Group A Streptococcal antigen in human throat swab specimens.

TEST PROCEDURE



RESULT INTERPRETATION



**RAPID TEST
STREP A**



www.prognosis-biotech.com

✉ exports@prognosis-biotech.com ☎ +30 2410 623922

📍 Farsalon 153 | 41335 Larissa, Greece

**PROGNOSIS
BIOTECH**