



Proficiency Test	Assigned value (mg/L)	Result (mg/L)	Z-score
FAPAS - Mustard in Durum Wheat - 27456A - December 2025 - January 2026	13.2	13.85	0.2
FAPAS - Mustard in Durum Wheat - 27456B - December 2025 - January 2026	22.2	24.12	0.4

Proficiency Test	Assignment	Result (mg/kg)	Satisfactory results / results agreeing with assignment (%)
FAPAS Celery and Mustard in Instant Soup Powder - 27434A - June - July 2025	DETECTED	DETECTED	75.0



FAPAS Celery and Mustard in Instant Soup Powder - 27434B - June - July 2025	DETECTED	DETECTED	100.0
---	----------	----------	-------

Proficiency Test	Assigned value (mg/kg) according to R-Biopharm Ridascreen Fast Senf/Mustard R6152	Range for $ z \leq 2$ according to R-Biopharm Ridascreen Fast Senf/Mustard R6152	Result (mg/kg)	Recovery (%)
FAPAS Celery & Mustard in Instant Soup Powder 27397A June - July 2024	43.4	21.7 - 65.1	44.0	101.4
FAPAS Celery & Mustard in Instant Soup Powder 27397B June - July 2024	NOT DETECTED	-	<1.0	-