

Ring Test Results

Proficiency Test	Assigned value (mg/kg) according to R-Biopharm RIDASCREEN FAST Lupine R6102	Range for $ z \leq 2$ according to R-Biopharm RIDASCREEN FAST Lupine R6102	Result (mg/kg)	Recovery (%)
FAPAS Lupin in Flour 27318A March – April 2022	19.4	9.71 – 29.1	18.55	95.6
FAPAS Lupin in Flour 27318B March – April 2022	NOT DETECTED	-	<1.0	-
FAPAS Lupin in Flour 27352A March – April 2023	23.8	11.9 – 35.7	31.5	132.4
FAPAS Lupin in Flour 27352B March – April 2023	NOT DETECTED	-	<1.0	-