

Proficiency Test	Assigned value (mg/kg) according to R-Biopharm Ridascreen Fast Soya R7102	Range for $ z \leq 2$ according to R-Biopharm Ridascreen Fast Soya R7102	Result (mg/kg)	Recovery (%)
FAPAS Soya in Wheat Flour 27356A - May – June 2023	NOT DETECTED (89% CONSENSUS)	-	<1.0	-
FAPAS Soya in Wheat Flour 27356B - May – June 2023	29.0	14.5 – 43.5	22.5	77.6
FAPAS Soya in Wheat Flour 27337A - September - November 2022	DETECTED (VALUE NOT SET)	-	9.18	-
FAPAS Soya in Wheat Flour 27337B - September - November 2022	20.9	10.4 – 31.3	17.38	83.2
FAPAS Soya in Wheat Flour 27323A - April – June 2022	DETECTED	-	DETECTED	-
FAPAS Soya in Wheat Flour 27323B - April – June 2022	15.2	7.58 – 22.8	21.0	138.2
FAPAS Soya in Wheat Flour 27295A - April – May 2021	10.4	5.19 – 15.6	11.49	110.5
FAPAS Soya in Wheat Flour 27295B - April – May 2021	26.0	13.0 – 39.0	26.32	101.2

Proficiency Test	Assignment	Result (mg/kg)	Satisfactory results / results agreeing with assignment (%)
FAPAS Allergens on Swab 27383E - January - March 2024	NOT DETECTED	NOT DETECTED	100.0
FAPAS Allergens on Swab 27383F - January - March 2024	DETECTED	DETECTED	63.0

Proficiency Test	Assigned value (mg/L)	Result (mg/L)	Z-score
BIPEA Food Allergens - Soya - 77H - Meat - 09 - 2277 - January 2024	DETECTED	DETECTED	-
BIPEA Food Allergens - Soya - 77H - Meat - 09 - 2377 - January 2024	DETECTED	DETECTED	-
BIPEA Food Allergens - Soya - 77H - Meat - 10 - 2477 - January 2024	NOT DETECTED	<0.2	-
BIPEA Food Allergens - Soya - 77H - Flour - 07 - 3177 - June 2023	NOT DETECTED	<1.0	-
BIPEA Food Allergens - Soya - 77H - Flour - 07 - 3277 - June 2023	64.0	42.8	-1.33
BIPEA Food Allergens - Soya - 77H - Flour - 08 - 3377 - June 2023	27.5	20.7	-0.88
BIPEA Food Allergens - Soya - 77H - Meat - 08 - 2277 - January 2023	DETECTED	DETECTED	-



BIPEA Food Allergens - Soya - 77H - Meat - 08 - 2377 - January 2023	NOT DETECTED	<0.2	-
BIPEA Food Allergens - Soya - 77H - Meat - 09 - 2477 - January 2023	DETECTED	DETECTED	-
BIPEA Food Allergens - Soya - 77H - Flour - 06 - 3177 - June 2022	15.7	13.2	-0.47
BIPEA Food Allergens - Soya - 77H - Flour - 06 - 3277 - June 2022	37.3	32.1	-0.62
BIPEA Food Allergens - Soya - 77H - Flour - 06 - 3377 - June 2022	NOT DETECTED	<1.0	-