



| Proficiency Test  | Assigned value (mg/L) | Result (mg/L) | Z-score |
|---|-----------------------|---------------|---------|
| FAPAS Gluten, Milk & Egg in Cooked Biscuit 27459 January - March 2026                                 | 19.9                  | 22.5          | 0.5     |
| BIPEA Food Allergens - Milk 77C - 20 - Infant Formula - 21 - 0777 - April 2025                        | NOT DETECTED          | <0.5          | -       |
| BIPEA Food Allergens - Milk 77C - 20 - Infant Formula - 21 - 0877 - April 2025                        | 26.6                  | DETECTED      | -       |
| BIPEA Food Allergens - Milk 77C - 20 - Infant Formula - 20 - 0977 - April 2025                        | 10.1                  | 10.1          | 0.0     |
| BIPEA Allergens in Industrial Process Water 77D - Industrial Process Water - 19 - 1077 - January 2025 | 9.9                   | 7.8           | -0.81   |
| BIPEA Allergens in Industrial Process Water 77D - Industrial Process Water - 18 - 1177 - January 2025 | NOT DETECTED          | <0.5          | -       |

|   |              |          |       |
|---|--------------|----------|-------|
| BIPEA Allergens in Industrial Process Water 77D - Industrial Process Water - 18 - 1277 - January 2025 | DETECTED     | DETECTED | -     |
| BIPEA Food Allergens - Milk 77C - 19 - Infant Formula - 20 - 0777 - November 2024                     | NOT DETECTED | <0.5     | -     |
| BIPEA Food Allergens - Milk 77C - 19 - Infant Formula - 20 - 0877 - November 2024                     | 2.0          | 1.6      | -1.0  |
| BIPEA Food Allergens - Milk 77C - 19 - Infant Formula - 19 - 0977 - November 2024                     | 15.7         | 13.4     | -0.55 |
| FAPAS Gluten, Egg and Milk in Cake Mix 27400A July - September 2024                                   | 14.4         | 14.7     | 0.1   |
| FAPAS Gluten, Egg and Milk in Cake Mix 27400B July - September 2024                                   | 30.4         | 31.2     | 0.1   |
| BIPEA Food Allergens - Milk 77C - 18 - Infant Formula - 19 - 0777 - April 2024                        | 20.5         | 25.8     | 1.03  |
| BIPEA Food Allergens - Milk 77C - 18 - Infant Formula - 19 - 0877 - April 2024                        | 33.6         | 35.8     | 0.26  |
| BIPEA Food Allergens - Milk 77C - 18 - Infant Formula - 18 - 0977 - April 2024                        | NOT DETECTED | <0.5     | -     |
| BIPEA Allergens in Wines 71 - 32 - White Wine - 11 - 0671 - February 2024                             | 0.51         | 0.36     | -1.30 |
| BIPEA Allergens in Industrial Process Water 77D - Industrial Process Water - 16 - 1077 - January 2024 | 30.8         | 26.9     | -0.32 |
| BIPEA Allergens in Industrial Process Water 77D - Industrial Process Water - 16 - 1177 - January 2024 | NOT DETECTED | <1.0     | -     |
| BIPEA Allergens in Industrial Process Water 77D - Industrial Process Water - 16 - 1277 - January 2024 | 34.0         | 27.6     | -0.43 |
| BIPEA Food Allergens - Milk 77C - 17 - 0280 Infant Formula - 18 - 0777 - November 2023                | 8.0          | 4.3      | -1.72 |
| BIPEA Food Allergens - Milk 77C - 17 - 0280 Infant Formula - 18 - 0877 - November 2023                | 4.1          | 2.2      | -1.19 |

|  |              |          |      |
|--|--------------|----------|------|
| BIPEA Food Allergens - Milk 77C - 17 - 0280 Infant Formula<br>- 17 - 0977 - November 2023                                    | NOT DETECTED | <1.0     | -    |
| BIPEA Food Allergens - Milk 77C - 15 - 0204 Infant Formula<br>- 16 - 0777 - November 2022                                    | 5.4          | 8.6      | 1.83 |
| BIPEA Food Allergens - Milk 77C - 15 - 0204 Infant Formula<br>- 17 - 0877 - November 2022                                    | 27.6         | 44.9     | 1.93 |
| BIPEA Food Allergens - Milk 77C - 15 - 0204 Infant Formula<br>- 15 - 0977 - November 2022                                    | NOT DETECTED | <2.5     | -    |
| BIPEA Allergens in Wines 71 - 28 - Rose Wine - 09 - 1271 -<br>October 2022   | 0.30         | <0.60    | -    |
| BIPEA Allergens in Industrial Process Water 77D - 0894<br>Industrial Process Water - 13 - 1077 - June 2022                   | NOT DETECTED | <2.5     | -    |
| BIPEA Allergens in Industrial Process Water 77D - 0894<br>Industrial Process Water - 13 - 1177 - June 2022                   | DETECTED     | DETECTED | -    |
| BIPEA Allergens in Industrial Process Water 77D - 13 2021-<br>2022 - 0894 Industrial Process Water - 13 -1277 - June<br>2022 | DETECTED     | DETECTED | -    |
| BIPEA Allergens in Wines 71 - 27 - 0902 - Red Wine - 08 -<br>0971- June 2022   | NOT DETECTED | <0.25    | -    |

| <b>Proficiency Test</b>   | <b>Assigned value<br/>(mg/kg) according<br/>to R-Biopharm<br/>Ridascreen Fast<br/>Casein R4612</b> | <b>Range for <math> z  \leq 2</math><br/>according to R-<br/>Biopharm<br/>Ridascreen Fast<br/>Casein R4612</b> | <b>Result<br/>(mg/kg)</b> | <b>Recovery<br/>(%)</b> |
|---|--|--|---------------------------|-------------------------|
| FAPAS Gluten, Milk & Egg in Cake<br>Mix 27438A July - August 2025     | 17.5   | 8.74 - 26.2  | 20.8                      | 118.9                   |
| FAPAS Gluten, Milk & Egg in Cake<br>Mix 27438B July - August 2025     | 26.5   | 13.3 - 39.8  | 39.6                      | 149.4                   |
| FAPAS Milk in Infant Breakfast<br>Cereal 27435A June - August<br>2025 | 14.0   | 7.01 - 21.0  | 13.7                      | 97.9                    |

|   |              |             |       |       |
|---|--------------|-------------|-------|-------|
| FAPAS Milk in Infant Breakfast Cereal 27435B June - August 2025           | NOT DETECTED | -           | <0.50 | -     |
| FAPAS Milk & Gluten in Infant Soya Formula 27422A February - March 2025   | 19.7         | 9.86 - 29.6 | 22.2  | 112.7 |
| FAPAS Milk & Gluten in Infant Soya Formula 27422B February - March 2025   | NOT DETECTED | -           | <0.50 | -     |
| FAPAS Gluten, Milk & Egg in Cooked Biscuit 27418 January - March 2025     | 17.9         | 8.93 - 26.8 | 17.8  | 99.4  |
| FAPAS Milk in Infant Breakfast Cereal - 27398A - June - August 2024       | NOT DETECTED | -           | <0.5  | -     |
| FAPAS Milk in Infant Breakfast Cereal - 27398B - June - August 2024       | 19.1         | 9.55 - 28.6 | 13.4  | 70.2  |
| FAPAS Milk and Gluten in Infant Soya Formula 27386A February - April 2024 | NOT DETECTED | -           | <0.5  | -     |
| FAPAS Milk and Gluten in Infant Soya Formula 27386B February - April 2024 | 13.3         | 6.65 - 20.0 | 17.7  | 133.1 |
| FAPAS Allergens in Cooked Biscuit 27382 January - March 2024              | 13.0         | 6.50 - 19.5 | 15.02 | 115.5 |
| FAPAS Allergens in Cake Mix - 27364A - July - September 2023              | 14.2         | 7.09 - 21.3 | 14.1  | 99.3  |
| FAPAS Allergens in Cake Mix - 27364B - July - September 2023              | 35.0         | 17.5 - 52.5 | 31.4  | 89.7  |
| FAPAS Milk in Infant Breakfast Cereal - 27362A - July - August 2023       | NOT DETECTED | -           | <1.0  | -     |

|   |              |             |          |       |
|---|--------------|-------------|----------|-------|
| FAPAS Milk in Infant Breakfast Cereal - 27362B - July - August 2023         | 18.0         | 9.00 - 27.0 | 16.1     | 89.4  |
| FAPAS Allergens in Cake Mix - 27331A - July - September 2022                | 12.4         | 6.18 - 18.5 | 10.75    | 86.7  |
| FAPAS Allergens in Cake Mix - 27331B - July - September 2022                | 5.63         | 2.82 - 8.45 | 4.76     | 84.5  |
| FAPAS Milk in Infant Breakfast Cereal - 27329A - June - August 2022         | 15.2         | 7.58 - 22.7 | 20.1     | 132.2 |
| FAPAS Milk in Infant Breakfast Cereal - 27329B - June - August 2022         | NOT DETECTED | -           | <1.0     | -     |
| FAPAS Gluten & Milk in Infant Soya Formula - 27316A - February - April 2022 | 31.0         | 15.5 - 46.5 | 23.04    | 74.3  |
| FAPAS Gluten & Milk in Infant Soya Formula - 27316B - February - April 2022 | NOT DETECTED | -           | <1.0     | -     |
| FAPAS Allergens in Cooked Biscuit - 27314 - January - March 2022            | DETECTED     | -           | DETECTED | -     |
| FAPAS Milk in Infant Breakfast Cereal - 27322A - September - October 2021   | NOT DETECTED | -           | <1.0     | -     |
| FAPAS Milk in Infant Breakfast Cereal - 27322B - September - October 2021   | 21.1         | 10.5 - 31.6 | 19.59    | 92.8  |
| FAPAS Milk in Infant Formula (hydrolyzed) - 27304A - August - October 2021  | NOT DETECTED | -           | <1.0     | -     |
| FAPAS Milk in Infant Formula (hydrolyzed) - 27304B - August - October 2021  | DETECTED     | -           | DETECTED | -     |



|  |      |             |       |       |
|--|------|-------------|-------|-------|
| FAPAS Allergens in Cake Mix - 27301A - July - September 2021     | 11.0 | 5.5 - 16.5  | 13.99 | 127.2 |
| FAPAS Allergens in Cake Mix - 27301B - July - September 2021     | 6.0  | 3.0 - 9.0   | 5.2   | 86.7  |
| FAPAS Allergens in Cooked Biscuit - 27286 - January - March 2021 | 30.1 | 15.0 - 45.1 | 28.7  | 95.3  |