

Proficiency Test	Assignment	Result (mg/kg)	Satisfactory results / results agreeing with assignment (%)
FAPAS Fish and Peanut in Sauce (Cooked) - 27412A - December 2024 - February 2025	DETECTED	DETECTED	56.0

Proficiency Test	Assigned value (mg/kg) according to Romer Labs Agraquant Fish (10002083)	Range for $ z \leq 2$ according to Romer Labs Agraquant Fish (10002083)	Result (mg/kg)	Recovery (%)
FAPAS - Fish and Peanut in Sauce (Cooked) - 27451A - December 2025 - February 2026	30.8	15.4 - 46.2	23.2	75.3
FAPAS - Fish and Peanut in Sauce (Cooked) - 27451B - December 2025 - February 2026	NOT DETECTED	-	<1.0	-
FAPAS Fish in Sauce (Cooked) - 27376A - November 2023 - January 2024	NOT DETECTED	-	<2.5	-
FAPAS Fish in Sauce (Cooked) - 27376B - November 2023 - January 2024	34.2	17.1 - 51.3	27.4	80.1