




Proficiency Test	Assigned value (mg/kg) of allergenic protein	Result (mg/kg)	Z-score
FAPAS Lupin in Flour 27465A March – May 2026	NOT DETECTED	<1.0	-
FAPAS Lupin in Flour 27465B March – May 2026	17.8	17.8	0.0
FAPAS Lupin in Flour 27388A February – April 2024	19.6	18.1	-0.3
FAPAS Lupin in Flour 27388B February – April 2024	NOT DETECTED	<1.0	-

Proficiency Test	Assigned value (mg/kg) according to R-Biopharm RIDASCREEN FAST Lupine R6102	Range for $z \leq 2$ according to R-Biopharm RIDASCREEN FAST Lupine R6102	Result (mg/kg)	Recovery (%)
FAPAS Lupin in Flour 27424A March – April 2025	11.4	5.68 – 17.0	12.12	106.3
FAPAS Lupin in Flour 27424B March – April 2025	NOT DETECTED	-	<1.0	-
FAPAS Lupin in Flour 27388A February – April 2024	20.9	10.4 – 31.3	18.1	86.6
FAPAS Lupin in Flour 27388B February – April 2024	NOT DETECTED	-	<1.0	-
FAPAS Lupin in Flour 27352A March – April 2023	23.8	11.9 – 35.7	31.5	132.4
FAPAS Lupin in Flour 27352B March – April 2023	NOT DETECTED	-	<1.0	-
FAPAS Lupin in Flour 27318A March – April 2022	19.4	9.71 – 29.1	18.55	95.6
FAPAS Lupin in Flour 27318B March – April 2022	NOT DETECTED	-	<1.0	-