



Proficiency Test	Assigned value (mg/L)	Result (mg/L)	Z-score
FAPAS - Mustard in Durum Wheat - 27456A - December 2025 - January 2026	51.5	52.9	0.1
FAPAS - Mustard in Durum Wheat - 27456B - December 2025 - January 2026	104	106	0.1

Proficiency Test	Assignment	Result (mg/kg)	Satisfactory results / results agreeing with assignment (%)
FAPAS Celery and Mustard in Instant Soup Powder - 27434A - June - July 2025	DETECTED	DETECTED	75.0

FAPAS Celery and Mustard in Instant Soup Powder - 27434B - June - July 2025	DETECTED	DETECTED	100.0
---	----------	----------	-------

Proficiency Test	Assigned value (mg/kg) according to R-Biopharm Ridascreen Fast Senf/Mustard R6152	Range for $z \leq 2$ according to R-Biopharm Ridascreen Fast Senf/Mustard R6152	Result (mg/kg)	Recovery (%)
FAPAS Celery & Mustard in Instant Soup Powder 27397A June - July 2024	43.4	21.7 - 65.1	41.5	95.6
FAPAS Celery & Mustard in Instant Soup Powder 27397B June - July 2024	NOT DETECTED	-	<1.0	-
FAPAS Celery & Mustard in Instant Soup Powder 27361A June - July 2023	44.1	22.0 - 66.1	46.1	104.5
FAPAS Celery & Mustard in Instant Soup Powder 27361B June - July 2023	NOT DETECTED	-	<1.0	-
FAPAS Celery & Mustard in Instant Soup Powder 27328A June - August 2022	NOT DETECTED	-	<1.0	-
FAPAS Celery & Mustard in Instant Soup Powder 27328B June - August 2022	DETECTED	-	DETECTED	-