

| Proficiency Test | Assigned value (mg/L or mg/kg) | Result (mg/L or mg/kg) | Z-score |
|---|--------------------------------|------------------------|---------|
| FAPAS Gluten & Sesame in Cumin 27423B February - April 2025 | 3.75 | 3.75 | 0.0 |

| Proficiency Test | Assigned value (mg/kg) according to R-Biopharm Ridascreen Fast Sesame R7202 | Range for $ z \leq 2$ according to R-Biopharm Ridascreen Fast Sesame R7202 | Result (mg/kg) | Recovery (%) |
|---|---|---|----------------|--------------|
| FAPAS Gluten & Sesame in Spice - Cumin 27464A February - April 2026 | 36.7 | 18.3 - 55.0 | 47.7 | 130.0 |
| FAPAS Gluten & Sesame in Spice - Cumin 27464B February - April 2026 | NOT DETECTED | - | <1.0 | - |

| Proficiency Test | Assignment | Result (mg/kg) | Satisfactory results / results agreeing with assignment (%) |
|---|--------------|----------------|---|
| FAPAS Gluten & Sesame in Cumin 27423A February - April 2025 | DETECTED | DETECTED | 100 |
| FAPAS Gluten & Sesame in Cumin 27423B February - April 2025 | DETECTED | DETECTED | 97 |
| FAPAS Sesame in Bread 27409B September - November 2024 | NOT DETECTED | <1.0 | 96 |
| FAPAS Sesame in Bread 27409A September - November 2024 | DETECTED | DETECTED | 96 |



| | | | |
|---|--------------|----------|-----|
| FAPAS Sesame in Bread 27373A October - November 2023 | NOT DETECTED | <1.0 | 100 |
| FAPAS Sesame in Bread 27373B October - November 2023 | DETECTED | DETECTED | 96 |
| FAPAS Gluten & Sesame in Cumin 27317A February - April 2022 | DETECTED | DETECTED | 100 |
| FAPAS Gluten & Sesame in Cumin 27317B February - April 2022 | NOT DETECTED | <1.0 | 82 |