



| Proficiency Test | Assigned value (mg/kg) | Result (mg/kg) | Z-score |
|---|------------------------|----------------|---------|
| BIPEA Food Allergens - Soya - 77H - Meat - 11 - 2277 - January 2026 | NOT DETECTED | <2.5 | - |
| BIPEA Food Allergens - Soya - 77H - Meat - 11 - 2377 - January 2026 | 6.40 | 10.8 | 1.96 |
| BIPEA Food Allergens - Soya - 77H - Meat - 11 - 2477 - January 2026 | 32.7 | 39.8 | 0.70 |
| FAPAS Soya in Wheat Flour 27446A - September - October 2025 | NOT DETECTED | <1.0 | - |
| FAPAS Soya in Wheat Flour 27446B - September - October 2025 | 30.3 | 24.7 | -0.7 |

| | | | |
|---|--------------|------|-------|
| BIPEA Food Allergens - Soya - 77H - Wheat Flour - 09 - 3177 - June 2025 | 50.4 | 44.6 | -0.61 |
| BIPEA Food Allergens - Soya - 77H - Wheat Flour - 09 - 3277 - June 2025 | NOT DETECTED | <2.5 | - |
| BIPEA Food Allergens - Soya - 77H - Wheat Flour - 10 - 3377 - June 2025 | 18.3 | 19.2 | 0.60 |
| BIPEA Food Allergens - Soya - 77H - Meat - 10 - 2277 - January 2025 | 58.4 | 82.4 | 1.14 |
| BIPEA Food Allergens - Soya - 77H - Meat - 10 - 2377 - January 2025 | NOT DETECTED | <2.5 | - |
| BIPEA Food Allergens - Soya - 77H - Meat - 10 - 2477 - January 2025 | 27.2 | 38.4 | 1.07 |
| BIPEA Food Allergens - Soya - 77H - Wheat Flour - 09 - 3377 - June 2024 | 39.1 | 54.1 | 1.32 |
| BIPEA Food Allergens - Soya - 77H - Wheat Flour - 08 - 3277 - June 2024 | NOT DETECTED | <2.5 | - |
| BIPEA Food Allergens - Soya - 77H - Wheat Flour - 08 - 3177 - June 2024 | 14.1 | 18.7 | 1.18 |
| FAPAS Soya in Wheat Flour 27392B - April - May 2024 | 21.9 | 24.5 | 0.5 |
| BIPEA Food Allergens - Soya - 77H - Flour - 07 - 3177 - June 2023 | NOT DETECTED | <2.5 | - |
| BIPEA Food Allergens - Soya - 77H - Flour - 07 - 3277 - June 2023 | 64.0 | 42.8 | -1.33 |
| BIPEA Food Allergens - Soya - 77H - Flour - 08 - 3377 - June 2023 | 27.5 | 20.7 | -0.88 |
| BIPEA Food Allergens - Soya - 77H - Flour - 06 - 3177 - June 2022 | 15.7 | 13.2 | -0.47 |

| | | | |
|---|--------------|------|-------|
| BIPEA Food Allergens - Soya - 77H - Flour - 06 - 3277 - June 2022 | 37.3 | 32.1 | -0.62 |
| BIPEA Food Allergens - Soya - 77H - Flour - 06 - 3377 - June 2022 | NOT DETECTED | <2.5 | - |

| Proficiency Test | Assigned value (mg/kg) according to R-Biopharm Ridascreen Fast Soya R7102 | Range for $z \leq 2$ according to R-Biopharm Ridascreen Fast Soya R7102 | Result (mg/kg) | Recovery (%) |
|--|--|---|-----------------------|---------------------|
| FAPAS Soya in Wheat Flour 27356A - May - June 2023 | NOT DETECTED (89% CONSENSUS) | - | <1.0 | - |
| FAPAS Soya in Wheat Flour 27356B - May - June 2023 | 29.0 | 14.5 - 43.5 | 22.5 | 77.6 |
| FAPAS Soya in Wheat Flour 27337A - September - November 2022 | DETECTED (VALUE NOT SET) | - | 9.18 | - |
| FAPAS Soya in Wheat Flour 27337B - September - November 2022 | 20.9 | 10.4 - 31.3 | 17.38 | 83.2 |
| FAPAS Soya in Wheat Flour 27323A - April - June 2022 | DETECTED | - | DETECTED | - |
| FAPAS Soya in Wheat Flour 27323B - April - June 2022 | 15.2 | 7.58 - 22.8 | 21.0 | 138.2 |
| FAPAS Soya in Wheat Flour 27295A - April - May 2021 | 10.4 | 5.19 - 15.6 | 11.49 | 110.5 |
| FAPAS Soya in Wheat Flour 27295B - April - May 2021 | 26.0 | 13.0 - 39.0 | 26.32 | 101.2 |

| Proficiency Test | Assignment | Result (mg/kg) | Satisfactory results / results agreeing with assignment (%) |
|---|-------------------|-----------------------|--|
| FAPAS Allergens on Swab 27419E - May - July 2025 | NOT DETECTED | <0.05 | 91.0 |
| FAPAS Allergens on Swab 27419F - May - July 2025 | DETECTED | DETECTED | 39.0 |
| FAPAS Soya in Wheat Flour 27392A - April - May 2024 | DETECTED | DETECTED | 100.0 |
| FAPAS Soya in Wheat Flour 27392B - April - May 2024 | DETECTED | DETECTED | 100.0 |
| FAPAS Allergens on Swab 27383E - January - March 2024 | NOT DETECTED | NOT DETECTED | 100.0 |
| FAPAS Allergens on Swab 27383F - January - March 2024 | DETECTED | DETECTED | 63.0 |
| BIPEA Food Allergens - Soya - 77H - Meat - 09 - 2277 - January 2024 | DETECTED | DETECTED | 100.0 |
| BIPEA Food Allergens - Soya - 77H - Meat - 09 - 2377 - January 2024 | DETECTED | DETECTED | 100.0 |
| BIPEA Food Allergens - Soya - 77H - Meat - 10 - 2477 - January 2024 | NOT DETECTED | <0.2 | 100.0 |
| BIPEA Food Allergens - Soya - 77H - Meat - 08 - 2277 - January 2023 | DETECTED | DETECTED | 100.0 |
| BIPEA Food Allergens - Soya - 77H - Meat - 08 - 2377 - January 2023 | NOT DETECTED | <0.2 | 100.0 |
| BIPEA Food Allergens - Soya - 77H - Meat - 09 - 2477 - January 2023 | DETECTED | DETECTED | 100.0 |