




Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score
BIPEA 31E - Mycotoxins - Feed - Soyabean Meal - 18 - 6031 - May 2026	16.51	13.80	-0.62
BIPEA 31A - Mycotoxins - Dried Fruits, Spices and Other Products - Lentils - 03 - 5931 - May 2026	3.66	4.00	0.31
BIPEA 31B - Mycotoxins - Cereals - Ground Barley - 06 - 0531 - February 2026	2.64	4.11	1.86
BIPEA 31A - Mycotoxins - Dried Fruits, Spices and Other Products - Coffee - 06 - 1731 - January 2026	6.26	6.00	-0.14



FAPAS Food Chemistry Proficiency Test - 04578 – Multi-mycotoxins in Durum Wheat - November 2025 – January 2026	3.80	3.68	-0.1
FAPAS Food Chemistry Proficiency Test - 04575 – Aflatoxins in Hemp (CBD) leaf - November 2025 – January 2026	5.70	4.42	-1.0
BIPEA 31A - Mycotoxins – Dried Fruits, Spices and Other Products – Thyme – 01- 6331 – November 2025	3.77	2.85	-0.81
FAPAS Food Chemistry Proficiency Test - 04559 – Multi - Mycotoxins in Pet Dog Food (Dry, Cereal Based) - July – August 2025	7.17	7.70	0.3
FAPAS Food Chemistry Proficiency Test - 04548 – Mycotoxins in Animal Feed (Cereal Based) - April – May 2025	5.57	3.38	-1.8
FAPAS Food Chemistry Proficiency Test - 04543 – Mycotoxins in Cereal Snack - February – April 2025	1.83	2.4	1.4
FAPAS Food Chemistry Proficiency Test - 04533 – Mycotoxin Contamination in Animal Feed (Cereal Based) - November 2024 – January 2025	3.59	3.6	0.0
FAPAS Food Chemistry Proficiency Test - 04546 - Aflatoxins in Hemp (CBD) Leaf - November 2024 – January 2025	3.66	5.02	1.7
BIPEA 31C - Mycotoxins – Baby Food – 24- 3931 – December 2024	3.91	5.01	0.94



FAPAS Food Chemistry Proficiency Test - 22231 - Multi - Mycotoxins in Wheat Flour - November 2024 – January 2025	3.92	4.01	0.1
FAPAS Food Chemistry Proficiency Test - 04528 - Mycotoxins in Animal Feed (Cereal Based) - September - November 2024	4.28	4.56	0.3
BIPEA 31B - Mycotoxins – Cereals – Wheat Flour -1- 6231 – October 2024	3.45	5.51	1.99
BIPEA 31C - Mycotoxins – Baby Food – 23 - 3931 – June 2024	7.50	8.16	0.31
BIPEA 31A - 188 - Mycotoxins – Dried Fruits, Spices and Other Products – Hazelnut Paste – 3 - 2031 – March 2024	5.03	6.91	1.25
BIPEA 31B – 7 - 188 - Mycotoxins – Cereals – Oats – 7 - 2931 – February 2024	1.21	1.62	1.12
BIPEA 31A – 5 - 187 - Mycotoxins – Dried Fruits, Spices and Other Products – Coffee – 5 - 1731 – January 2024	7.89	6.38	-0.68
FAPAS Food Chemistry Proficiency Test - 04500 - Multi - Mycotoxins in Durum Wheat - December 2023 – January 2024	2.31	2.23	-0.2
FAPAS Food Chemistry Proficiency Test - 22215 - Multi - Mycotoxins in Wheat Flour - November 2023 – January 2024	1.91	2.27	0.9
FAPAS Food Chemistry Proficiency Test - 04496 - Aflatoxins in Rice – November 2023 – January 2024	4.13	3.19	-1.0



BIPEA 31E - 186 - Mycotoxins - Feed – Corn Silage – 15 - 5831 - November 2023	5.48	4.64	-0.52
BIPEA 31A - 186 - Mycotoxins - Dried Fruits, Spices and Other Products – Ground Ginger – 02 - 1231 - November 2023	7.97	8.00	0.01
BIPEA 31B – 186 – Mycotoxin Contamination in Cereals – Common Wheat – 3 - 0631 - October 2023	3.30	2.63	-0.68
FAPAS Food Chemistry Proficiency Test - 04481 - Mycotoxins in Pet Dog Food - July – August 2023	9.88	5.66	-1.9
BIPEA 31B – 185 – Mycotoxin Contamination in Cereals – Rice - 06-4331 - June 2023	4.47	3.14	-0.99
BIPEA 31A - 184 Mycotoxins - Dried Fruits, Spices and Other Products – Pistachio Paste – 08 - 1431 - May 2023	10.05	13.45	1.23
A.I.A. Ring Test Lotto - RT B1 300523 - Aflatoxin B1 in Corn - May 2023	3.51	2.76	-0.57
	7.92	12.46	1.57
	15.08	5.86	-1.80
	27.21	26.00	-0.16
BIPEA 31B – 184 – Mycotoxin Contamination in Cereals – Rusk (Crispbread) 02-5431 - April 2023	2.65	3.63	1.23



BIPEA 31E - 3 - Mycotoxin Contamination in Animal feed - 13 - 3031 - January 2023		3.16	3.65	0.52
BIPEA 31B - 183 - Mycotoxin Contamination in Ground Barley - 05 - 0531 - December 2022		6.42	6.62	0.11
BIPEA 31B - 182 - Mycotoxin Contamination in Ground Corn - 18 - 0731 - December 2022		5.7	2.9	-1.65
BIPEA 31B - 180 - Mycotoxin Contamination in Rice - 05 - 4331 - June 2022		2.64	3.65	1.28
BIPEA 31E - 2 - Mycotoxin Contamination in Corn Silage - 01 - 5831 - May 2022		2.16	1.12	-1.6
FAPAS Food Chemistry Proficiency Test - 04430 - Mycotoxin Contamination in Animal Feed - December 2021 - January 2022		19.1	19.5	0.5
QUALINOVA ROUND : AL 21008 - Aflatoxins in dried fruit - April 2021		3.95	3.29	-0.76
DRRR - RVEP 190663 - CORN AFLATOXIN 1	Sample 1	8.6	9.4	0.53
	Sample 2	15.2	15.1	-0.03
FAPAS Food Chemistry Proficiency Test - 04366 - Multi - Mycotoxins in Maize - August - September 2019		1.88	1.84	-0.1



FAPAS Food Chemistry Proficiency Test - 04359 - Multi - Mycotoxins in Maize - April - May 2019	9.22	9.58	0.2
A.I.A. Ring Test Lotto - RT B1 260319 - Aflatoxin B1 in Corn - April 2019	<limit of detection	<0.5	0.117
	7.60	8.94	
	14.10	15.69	
	29.40	28.11	
FAPAS Food Chemistry Proficiency Test - 04354 - Multi - Mycotoxins in Maize - February - March 2019	4.85	4.58	-0.3
FAPAS Food Chemistry Proficiency Test - 04330 - Multi - Mycotoxins in Maize - February - March 2018	3.83	3.39	-0.5
FAPAS Food Chemistry Proficiency Test - 04319 - Multi - Mycotoxins in Maize - August - September 2017	4.69	5.12	0.4