



**BIO-SHIELD**

**BIO-SHIELD CHLORAMPHENICOL**

Cat.Number B8048/B8096



Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score
BIPEA 88C – 6 Veterinary Drugs in Meat Products 6 - 0388 – Meat - November 2023	21.0	11.3	-1.54
BIPEA 88B – 10 Veterinary Drugs in Honey 10 - 0288 - November 2023	36.5	27.1	-0.86
BIPEA 88C – 5 Veterinary Drugs in Meat Products 5-0388 - Meat - April 2023	20.2	25.5	0.88
BIPEA 88B – 9 Veterinary Drugs in Honey 9 - 0288 - April 2023	25.2	25.6	0.05



**PROGNOSIS**  
**BIOTECH**



## BIO-SHIELD

FAPAS 02488 – Phenicol in Prawns - November - December 2022	0.694	0.437	-1.7
BIPEA 88C – 4 Veterinary Drugs in Meat Products 4 - 0388 - Meat - November 2022	6.7	7.7	0.50
BIPEA 88B – 8 Veterinary Drugs in Honey 8 - 0288 - November 2022	15.6	9.9	-1.21