




BIO-SHIELD T-2/HT-2 8

Cat.Number B5648/B5696



Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score
BIPEA 31E - Mycotoxins - Feed - Soyabean Meal - 18 - 6031 - May 2026	104	110	0.22
FAPAS Food Chemistry Proficiency Test - 04578 - Multi-mycotoxins in Durum Wheat - November 2025 - January 2026	118	110	-0.30
BIPEA 31B - Mycotoxins - Cereals - Ground Corn - 19 - 0731 - December 2025	139	115	-0.66
FAPAS Food Chemistry Proficiency Test - 04559 - Multi - Mycotoxins in Pet Dog Food (Dry, Cereal Based) - July - August 2025	47.7	59.1	1.1



BIO-SHIELD

BIPEA 31B - Mycotoxins - Cereals - Ground Rye - 4 - 5131 - June 2025	35	39	0.38
BIPEA 31E - Mycotoxins - Feed - Animal Feed - 15 - 3031 - January 2025	48	62	1.04
BIPEA 31C - Mycotoxins - Baby Food - 24-3931 - December 2024	76	79	0.15
BIPEA 31B - Mycotoxins - Cereals - Maize Flour - 4 - 2231 - December 2024	102	99.0	-0.11
BIPEA 31B - Mycotoxins - Cereals - Wheat Flour - 1 - 6231 - October 2024	115	83	-1.05
FAPAS Food Chemistry Proficiency Test - 04520 - Multi - Mycotoxins in Maize Flour - July - September 2024	96.8	109	0.6
BIPEA 31B - Mycotoxins - Cereals - Triticale - 3 - 1931 - June 2024	79	70	-0.42
FAPAS Food Chemistry Proficiency Test - 04510 - Multi - Mycotoxins in Maize Flour - April - May 2024	111	97.23	-0.6
BIPEA 31B - Cereals - Bran - 4 - 3631 - April 2024	102	140	1.41
BIPEA 31B - 7 - 188 - Mycotoxins - Cereals - Oats - 7 - 2931 - February 2024	223	209	-0.24
BIPEA 31E - Mycotoxins - Feed - Animal Feed - 14 - 5 - January 2024	94	85	-0.36





BIO-SHIELD

FAPAS Food Chemistry Proficiency Test - 04500 - Multi - Mycotoxins in Durum Wheat - December 2023 - January 2024	80.9	94.0	0.7
BIPEA 31B - Mycotoxins - Cereals - Ground Corn - 6 - 187 - December 2023	210	244	0.63
BIPEA 31B - 186 - Mycotoxin Contamination in Cereals - Common Wheat - 3 - 0631 - October 2023	86.0	64.0	-0.96
FAPAS Food Chemistry Proficiency Test 04482 - Multi-Mycotoxins in Maize Flour - July - September 2023	104	103.2	0.0
FAPAS Food Chemistry Proficiency Test 04481 - Mycotoxins in Pet Dog Food - July - August 2023	40.2	39.1	-0.1
BIPEA 31c 20 Mycotoxins - Baby Food - 21 - 3931 - June 2023	125	129	0.12
BIPEA 31B - 185 Mycotoxin Contamination in Cereals - Rice - 06-4331 - June 2023	139	103	-0.99
FAPAS Food Chemistry Proficiency Test 04472 - Multi-Mycotoxins in Maize Flour - April - June 2023	169	166.7	-0.1
BIPEA 31B-183 Mycotoxin Contamination in Ground barley - 05-0531 - December 2022	84.0	83.0	-0.04
BIPEA 31E-3 Mycotoxin Contamination in Animal feed - 13-3031 - January 2023	113	137	0.8





FAPAS Food Chemistry Proficiency Test 22199 - Multi-Mycotoxins in Wheat Flour - November 2022-January 2023	93.7	72.54	-1.0
BIPEA 31B-182 Mycotoxin Contamination in Ground corn - 18-0731 - December 2022	116	103	-0.43
FAPAS Food Chemistry Proficiency Test 04448 - Multi-Mycotoxins in Maize Flour - July - September 2022	92.6	123.42	1.5
BIPEA 31B-180 Mycotoxin Contamination in Rice - 05-4331 - June 2022	99	66	-1.25
FAPAS Food Chemistry Proficiency Test 04439 - Multi-Mycotoxins in Maize Flour - April - May 2022	189	190.85	0.0
FAPAS Food Chemistry Proficiency Test 22187 - Multi-Mycotoxins in Wheat Flour - November 2021 - January 2022	54	61.1	0.6
FAPAS Food Chemistry Proficiency Test 04418 - Multi - Mycotoxins in Maize Flour - July - September 2021	111	122.64	0.5
FAPAS Food Chemistry Proficiency Test 04411 - Multi-Mycotoxins in Maize Flour - April - June 2021	174	169.47	-0.1
FAPAS Food Chemistry Proficiency Test 22177 - Multi - mycotoxins in Wheat Flour - November 2020 - January 2021	53.5	65.24	1.0
FAPAS Food Chemistry Proficiency Test 22175 - Deoxynivalenol, Zearalenone, T-2	523	511.96	-0.1



BIO-SHIELD

and HT-2 in Animal Feed - September - October 2020			
FAPAS Food Chemistry Proficiency Test 04384 - Multi Mycotoxins in Maize Flour - April - May 2020	330	247.37	-1.3

