



Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score
BIPEA 31A - Mycotoxins – Dried Fruits, Spices and Other Products – Almond Powder – 05 - 1131 – March 2026	8.20	7.04	-0.50
BIPEA 31B - Mycotoxins – Cereals – Ground Barley – 06- 0531 – February 2026	8.70	7.00	-0.70
BIPEA 31A - Mycotoxins – Dried Fruits, Spices and Other Products – Coffee – 06-1731 – January 2026	21.33	18.0	-0.60

BIO-SHIELD

FAPAS Food Chemistry Proficiency Test - 04575 - Aflatoxins in Hemp (CBD) leaf - November 2025 - January 2026	13.5	12.3	-0.4
BIPEA 31A - Mycotoxins - Dried Fruits, Spices and Other Products - Thyme - 01-6331 - November 2025	12.53	5.46	-1.77
BIPEA 31B - Mycotoxins - Cereals - Wheat Flour - 2 - 6231 - October 2025	3.91	3.82	-0.08
FAPAS Food Chemistry Proficiency Test - 04551 - Mycotoxin Contamination in Infant Food - April - June 2025	0.36	<0.5	-
FAPAS Food Chemistry Proficiency Test - 04543 - Mycotoxins in Cereal Snack - February - April 2025	4.19	4.07	-0.1
FAPAS Food Chemistry Proficiency Test - 04546 - Aflatoxins in Hemp (CBD) Leaf - November 2024 - January 2025	10.6	9.44	-0.5
BIPEA 31C - Mycotoxins - Baby Food -24-3931 - December 2024	14.50	16.50	0.52
FAPAS Food Chemistry Proficiency Test - 04528 - Mycotoxins in Animal Feed (Cereal Based) - September - November 2024	10.2	11.7	0.6
BIPEA 31C - Mycotoxins - Baby Food -23-3931 - June 2024	39.63	24.80	-1.46
BIPEA 31A - 188 - Mycotoxins - Dried Fruits, Spices and Other Products - Hazelnut Paste - 3 - 2031 - March 2024	19.51	20.60	0.21

BIO-SHIELD

BIPEA 31A - 186 Mycotoxins - Dried Fruits, Spices and Other Products - Ground Ginger - 02 - 1231 - November 2023			20.96	21.20	0.04
BIPEA 31A - Mycotoxins - Dried Fruits, Spices and Other Products - Peanut Paste - 6 - 185 - September 2023			8.69	9.57	0.36
BIPEA 31c 20 Mycotoxins - Baby Food - 21 - 3931 - June 2023			14.69	17.7	0.77
BIPEA 31A - 184 Mycotoxins Dried Fruits, Spices and Other Products - Pistachio Paste - 08-1431 May 2023			21.16	17.10	-0.73
QUALINOVA FOOD Mycotoxins in Feed - AL 23003 - April 2023	Sample 1		<1.0	0.989	-
	Sample 2		6.42	6.31	-0.07
BIPEA 31E-3 Mycotoxin Contamination in Animal feed 13 - 3031 January 2023			11.32	9.57	-0.57
BIPEA 31B-180 Mycotoxin Contamination in Rice 05 - 4331 June 2022			9.62	9.37	-0.09
BIPEA 31E-2 Mycotoxin Contamination in Corn Silage 01 -5831 May 2022			6.48	4.29	-1.17
FAPAS Food Chemistry Proficiency Test 04367 Mycotoxins in Dried Figs August - October 2019			10.8	8.1	-1.1
DRRR RVEP 190663 CORN AFLATOXIN 1	Sample 1		9.0	10.1	0.54
	Sample 2		17.9	16.6	-0.12



BIO-SHIELD

FAPAS Food Chemistry Proficiency Test 04364 Mycotoxins in Paprika June - August 2019	10.2	7.9	-1.0
FAPAS Food Chemistry Proficiency Test 04357 Aflatoxins in Peanut March - April 2019	12.40	13.86	0.5
FAPAS Food Chemistry Proficiency Test 04354 Multi Mycotoxins in Maize February - March 2019	4.95	4.78	-0.2
FAPAS Food Chemistry Proficiency Test 04330 Multi Mycotoxins in Maize February - March 2018	4.03	3.63	-0.4
BIPEA 31b 154 Mycotoxins - Cereals - Bran - 02-3631 - April 2017	28.0	21.5	-0.9