

### Ring Test Results

Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score	
BIPEA 31b 154 Mycotoxins - Cereals - Bran - 02-3631 - April 2017	28.0	21.5	-0.9	
FAPAS Food Chemistry Proficiency Test 04330 Multi Mycotoxins in Maize February - March 2018	4.03	3.63	-0.4	
FAPAS Food Chemistry Proficiency Test 04354 Multi Mycotoxins in Maize February - March 2019	4.95	4.78	-0.2	
FAPAS Food Chemistry Proficiency Test 04357 Aflatoxins in Peanut March - April 2019	12.40	13.86	0.5	
FAPAS Food Chemistry Proficiency Test 04364 Mycotoxins in Paprika June - August 2019	10.2	7.9	-1.0	
DRRR RVEP 190663 CORN AFLATOXIN 1	Sample 1	9.0	10.1	0.54
	Sample 2	17.9	16.6	-0.12
FAPAS Food Chemistry Proficiency Test 04367 Mycotoxins in Dried Figs August-October 2019	10.8	8.1	-1.1	
BIPEA 31E-2 Mycotoxin Contamination in Corn Silage 01-5831 May 2022	6.48	4.29	-1.17	



# BIO-SHIELD

BIPEA 31B-180 Mycotoxin Contamination in Rice 05-4331 June 2022	9.62	9.37	-0.09
--------------------------------------------------------------------	------	------	-------

