




Proficiency Test	Assigned value (mg/L)	Result (mg/L)	Z-score
BIPEA - Food Allergens - Nuts 77G – 19 - May 2025 - Cake - 13-2877	532.9	705.5	0.99
BIPEA - Food Allergens - Nuts 77G – 19 - May 2025 - Cake - 11-2977	318.6	450.8	1.13
BIPEA - Food Allergens - Nuts 77G – 19 - May 2025 - Cake - 15-3077	NOT DETECTED	<2.0	-



Proficiency Test	Assignment	Result (mg/kg)	Satisfactory results / results agreeing with assignment (%)
FAPAS - Hazelnut in Chocolate-27433A - May – July 2025	DETECTED	DETECTED	100.0
FAPAS - Hazelnut in Chocolate-27433B - May – July 2025	NOT DETECTED	<2.0	95.0
FAPAS - Hazelnut and Peanut in Cooked Biscuit -27425 - March – April 2025	DETECTED	DETECTED	95.0
FAPAS - Nuts In Chocolate - 27414 - November 2024 – January 2025	NOT DETECTED	NOT DETECTED	91.0
FAPAS – Hazelnut and Peanut in Cooked Biscuit - 27389 - March – May 2024	DETECTED	DETECTED	96.0
FAPAS - Nuts In Chocolate - 27378 - November 2023 – January 2024	NOT DETECTED	<2.0	100.0