

Proficiency Test	Assigned value (mg/L)	Result (mg/L)	Z-score
FAPAS Peanut as the Allergenic Protein in Chocolate 27437A - July – August 2025	3.76	4.0	0.3
FAPAS Peanut as the Allergenic Protein in Chocolate 27437B - July – August 2025	NOT DETECTED	<0.5	-

Proficiency Test	Assigned value (mg/kg) according to R-Biopharm Ridascreen Fast Peanut R6202	Range for $ z \leq 2$ according to R-Biopharm Ridascreen Fast Peanut R6202	Result (mg/kg)	Recovery (%)
FAPAS - Hazelnut and Peanut in Cooked Biscuit - 27425 - March – April 2025	34.5	17.3 – 51.8	22.6	65.5
FAPAS - Hazelnut and Peanut in Cooked Biscuit - 27389 - March – May 2024	7.29	3.64 – 10.9	9.8	134.4

Proficiency Test	Assignment	Result (mg/kg)	Satisfactory results / results agreeing with assignment (%)
FAPAS Allergens on Swab 27419A - May – July 2025	DETECTED	DETECTED	100.0
FAPAS Allergens on Swab 27419B - May – July 2025	NOT DETECTED	<0.5	92.0
FAPAS - Peanut in Chocolate – 27420A - February – March 2025	NOT DETECTED	<4.0	97.0



FAPAS - Peanut in Chocolate - 27420B - February - March 2025	DETECTED	DETECTED	97.0
FAPAS - Nuts In Chocolate - 27414 - November 2024 - January 2025	NOT DETECTED	<4.0	97.0
FAPAS - Hazelnut and Peanut in Cooked Biscuit -27389 March - May 2024	DETECTED	DETECTED	94.0
FAPAS Allergens on Swab 27383A - January - March 2024	DETECTED	DETECTED	69.0
FAPAS Allergens on Swab 27383B - January - March 2024	NOT DETECTED	<4.0	94.0