




Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score
UK Grain Testing Network - Mycotoxins Ring Check Sample [HY2026-R01-Tx01] [Wheat Flour] – July 2026	2.2	<35.0	-
UK Grain Testing Network - Mycotoxins Ring Check Sample [HY2026-R01-Tx02] [Wheat Flour] – July 2026	116.63	107.8	-0.2
UK Grain Testing Network - Mycotoxins Ring Check Sample [HY2026-R01-Tx03] [Wheat Flour] – July 2026	22.5	24.47	0.2
BIPEA 31B – Mycotoxins - Cereals – Pasta - 1 - 6431 – April 2026	61	65	0.24

FAPAS Food Chemistry Proficiency Test - 04578 - Multi-mycotoxins in Durum Wheat - November 2025 - January 2026	116	114	-0.1
UK Grain Testing Network - Mycotoxins Ring Check Sample [HY2025-R07-Tx01] [Wheat] - January 2026	85.86	69.3	-0.8
UK Grain Testing Network - Mycotoxins Ring Check Sample [HY2025-R07-Tx02] [Wheat] - January 2026	44.46	43	-0.1
UK Grain Testing Network - Mycotoxins Ring Check Sample [HY2025-R07-Tx03] [Wheat] - January 2026	1.8	<35.0	-
UK Grain Testing Network - Mycotoxins Ring Check Sample [HY2025-R04-Tx01] [Wheat] - December 2025	1.88	<35.0	-
UK Grain Testing Network - Mycotoxins Ring Check Sample [HY2025-R04-Tx02] [Wheat] - December 2025	48.2	39.0	-1.5
UK Grain Testing Network - Mycotoxins Ring Check Sample [HY2025-R04-Tx03] [Wheat] - December 2025	120.47	113.4	-0.3
BIPEA 31B - Mycotoxins - Cereals - Wheat Flour - 2 - 6231 - October 2025	25	26	0.13
UK NIR Grain Network Monthly Check Sample Results - HY2025-R01 - Wheat Flour (Spiked) - August 2025	19.2	22.0	0.3
UK NIR Grain Network Monthly Check Sample Results - HY2025-R02 - Wheat Flour (Spiked) - August 2025	108	110	-0.6

UK NIR Grain Network Monthly Check Sample Results - HY2025-R03 - Wheat Flour (Spiked) - August 2025	5	9	1.3
---	---	---	-----