

Proficiency Test	Assigned value (mg/L or mg/kg)	Result (mg/L or mg/kg)	Z-score
FAPAS Gluten & Sesame in Cumin 27423B February - April 2025	3.75	3.93	0.2

Proficiency Test	Assigned value (mg/kg) according to R-Biopharm Ridascreen Fast Sesame R7202	Range for $ z \leq 2$ according to R-Biopharm Ridascreen Fast Sesame R7202	Result (mg/kg)	Recovery (%)
FAPAS Gluten & Sesame in Spice - Cumin 27464A February - April 2026	36.7	18.3 - 55.0	41.4	112.8
FAPAS Gluten & Sesame in Spice - Cumin 27464B February - April 2026	NOT DETECTED	-	<2.0	-

Proficiency Test	Assignment	Result (mg/kg)	Satisfactory results / results agreeing with assignment (%)
FAPAS Gluten & Sesame in Cumin 27423A February - April 2025	DETECTED	DETECTED	100
FAPAS Gluten & Sesame in Cumin 27423B February - April 2025	DETECTED	DETECTED	97
FAPAS Sesame in Bread 27409B September - November 2024	NOT DETECTED	<4.0	96
FAPAS Sesame in Bread 27409A September - November 2024	DETECTED	DETECTED	96