



| Proficiency Test | Assigned value (mg/kg) | Result (mg/kg) | Z-score |
|---|------------------------|----------------|---------|
| BIPEA Food Allergens - Soya - 77H - Meat - 11 - 2277 - January 2026 | NOT DETECTED | <2.5 | - |
| BIPEA Food Allergens - Soya - 77H - Meat - 11 - 2377 - January 2026 | 6.40 | DETECTED | - |
| BIPEA Food Allergens - Soya - 77H - Meat - 11 - 2477 - January 2026 | 32.7 | 44.2 | 1.13 |
| FAPAS Soya in Wheat Flour 27446A - September - October 2025 | NOT DETECTED | <1.5 | - |
| FAPAS Soya in Wheat Flour 27446B - September - October 2025 | 30.3 | 26.9 | -0.5 |



| | | | |
|---|--------------|------|------|
| BIPEA Food Allergens - Soya - 77H - Wheat Flour - 09 - 3177 - June 2025 | 50.4 | 60.6 | 1.07 |
| BIPEA Food Allergens - Soya - 77H - Wheat Flour - 09 - 3277 - June 2025 | NOT DETECTED | <1.5 | - |
| BIPEA Food Allergens - Soya - 77H - Wheat Flour - 10 - 3377 - June 2025 | 18.3 | 18.8 | 0.33 |
| FAPAS Soya in Wheat Flour 27392B - April - May 2024 | 21.9 | 25.0 | 0.6 |

| Proficiency Test | Assignment | Result (mg/kg) | Satisfactory results / results agreeing with assignment (%) |
|---|--------------|----------------|---|
| FAPAS Allergens on Swab 27419E - May - July 2025 | NOT DETECTED | <0.1 | 91.0 |
| FAPAS Allergens on Swab 27419F - May - July 2025 | DETECTED | DETECTED | 39.0 |
| BIPEA Food Allergens - Soya - 77H - Meat - 10 - 2277 - January 2025 | DETECTED | DETECTED | 92.0 |
| BIPEA Food Allergens - Soya - 77H - Meat - 10 - 2377 - January 2025 | NOT DETECTED | NOT DETECTED | 92.0 |
| BIPEA Food Allergens - Soya - 77H - Meat - 10 - 2477 - January 2025 | DETECTED | DETECTED | 92.0 |
| FAPAS Soya in Wheat Flour 27392A - April - May 2024 | DETECTED | DETECTED | 100.0 |
| FAPAS Soya in Wheat Flour 27392B - April - May 2024 | DETECTED | DETECTED | 100.0 |
| FAPAS Allergens on Swab 27383E - January - March 2024 | NOT DETECTED | NOT DETECTED | 100.0 |



| | | | |
|---|----------|----------|------|
| FAPAS Allergens on Swab 27383F - January - March 2024 | DETECTED | DETECTED | 63.0 |
|---|----------|----------|------|