




Test	Assinged value (µg/kg)	Result (µg/kg)	Z-score
BIPEA 31B - Mycotoxins – Cereals – Ground Corn – 19 - 0731 – December 2025	259	300	0.62
FAPAS Food Chemistry Proficiency Test - 04549 – Multi - Mycotoxins in Maize Flour - April – May 2025	544	631	0.9
FAPAS Food Chemistry Proficiency Test - 04510 - Multi - Mycotoxins in Maize Flour - April – May 2024	788	803	0.1
BIPEA 31E - Mycotoxins – Feed – Animal Feed - 14 - 5 – January 2024	150	180	0.77

FAPAS Food Chemistry Proficiency Test - 04500 - Multi - Mycotoxins in Durum Wheat - December 2023 - January 2024	467	546	0.9
BIPEA 31B - Mycotoxins - Cereals - Ground Corn - 6 - 187 - December 2023	434	410	-0.22
FAPAS Food Chemistry Proficiency Test - 04482 - Multi - Mycotoxins in Maize Flour - July - September 2023	640	716	0.7
FAPAS Food Chemistry Proficiency Test - 04472 - Multi - Mycotoxins in Maize Flour - April - June 2023	798	785	-0.1
FAPAS Food Chemistry Proficiency Test - 04461 - Mycotoxin Contamination in Animal Feed (Cereal Based) - December 2022 - January 2023	347	323.67	-0.4
FAPAS Food Chemistry Proficiency Test - 04448 - Multi - Mycotoxins in Maize Flour - July - September 2022	569	563.73	-0.1
FAPAS Food Chemistry Proficiency Test - 04439 - Multi - Mycotoxins in Maize Flour - April - May 2022	983	920.56	-0.4
FAPAS Food Chemistry Proficiency Test - 04418 - Multi - Mycotoxins in Maize Flour - July - September 2021	720	705.87	-0.1
FAPAS Food Chemistry Proficiency Test - 04411 - Multi - Mycotoxins in Maize Flour - April - June 2021	962	951.18	-0.1



FAPAS Food Chemistry Proficiency Test - 04384 – Multi - Mycotoxins in Maize - April - May 2020	752	617.57	-1.1
FAPAS Food Chemistry Proficiency Test - 04366 – Multi - Mycotoxins in Maize - August – September 2019	650	629.2	-0.2
FAPAS Food Chemistry Proficiency Test - 04359 – Multi - Mycotoxins in Maize - April - May 2019	554	519.96	-0.4