

### Ring Test Results

Test	Assigned value (µg/kg)	Result (µg/kg)	Z-score
FAPAS Food Chemistry Proficiency Test 04359 Multi-Mycotoxins in Maize April - May 2019	554	519.96	-0.4
FAPAS Food Chemistry Proficiency Test 04366 Multi-Mycotoxins in Maize August - September 2019	650	629.2	-0.2
FAPAS Food Chemistry Proficiency Test 04384 Multi-Mycotoxins in Maize April - May 2020	752	617.57	-1.1
FAPAS Food Chemistry Proficiency Test 04411 Multi-Mycotoxins in Maize Flour April - June 2021	962	951.18	-0.1
FAPAS Food Chemistry Proficiency Test 04418 Multi - Mycotoxins in Maize Flour July - September 2021	720	705.87	-0.1
FAPAS Food Chemistry Proficiency Test 04439 Multi-Mycotoxins in Maize Flour April - May 2022	983	920.56	-0.4
FAPAS Food Chemistry Proficiency Test 04448 Multi-Mycotoxins in Maize Flour July - September 2022	569	563.73	-0.1