




Test	Assinged value (mg/kg)	Result (mg/kg)	Z-score
FAPAS Food Chemistry Proficiency Test 27467 Histamine in Fish (heat treated can/pouch) - March - May 2026	174	188.1	1.1
BIPEA Food Allergens - Histamine - 77E - Food Allergens - Histamine in Fish 21-1377 - March 2026	12.5	12.0	-0.2
BIPEA Food Allergens - Histamine - 77E - Food Allergens - Histamine in Fish 22-1477 - March 2026	NOT DETECTED	<8.0	-

BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 23-1577 – March 2026	34.7	33.2	-0.25
FAPAS Food Chemistry Proficiency Test 27441 Histamine in Fish (heat treated can/pouch) - August - September 2025	23.2	21.4	-0.8
Lab Conformance Dosage De l'histamine et de l'ABVT dans les produits de la pêche – AR-25 A - Niveau 1 – Juin 2025	22.5	19.4	-1.0
Lab Conformance Dosage De l'histamine et de l'ABVT dans les produits de la pêche – AR-25 A - Niveau 2 – Juin 2025	75.1	80.3	1.2
FAPAS Food Chemistry Proficiency Test 27426 Histamine in Canned Fish - March – May 2025	81.8	82.0	0.0
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 21-1577 – March 2025	NOT DETECTED	<8.0	-
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 20-1477 – March 2025	22.3	20.4	-0.38
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 19-1377 – March 2025	51.4	55.9	0.42
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 20-1577 - October 2024	NOT DETECTED	<8.0	-

BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 19-1477 - October 2024	45.2	35.7	-0.67
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 18-1377 - October 2024	25.8	20.1	-0.67
FAPAS Food Chemistry Proficiency Test 27403 Histamine in Canned Fish - August – September 2024	36.0	33.6	-0.7
FAPAS Food Chemistry Proficiency Test 27390 Histamine in Canned Fish - March – May 2024	273	262	-0.6
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 17-1377 - March 2024	94.5	91.8	-0.19
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 18-1477 - March 2024	NOT DETECTED	<8.0	-
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 19-1577 - March 2024	28.7	30.6	0.17
Lab Conformite Dosage De l'histamine et de l'ABVT dans les produits de la peche – AR-23 B - Niveau 1 - Octobre 2023	33.0	33.2	0.1
Lab Conformite Dosage De l'histamine et de l'ABVT dans les produits de la peche – AR-23 B - Niveau 2 - Octobre 2023	81.0	86.0	1.2

BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 16-1377 - October 2023	NOT DETECTED	<8.0	-
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 17-1477 - October 2023	61.5	61.1	-0.03
BIPEA Food Allergens - Histamine - 77E - Food Allergens – Histamine in Fish 18-1577 - October 2023	21.4	22.5	0.25
FAPAS Food Chemistry Proficiency Test 27367 Histamine in Canned Fish - August – September 2023	27.6	27.1	-0.2
Lab Conformance Dosage De l'histamine et de l'ABVT dans les produits de la pêche – Niveau 1 - Juin 2023	22.7	19.5	-0.8
Lab Conformance Dosage De l'histamine et de l'ABVT dans les produits de la pêche – Niveau 2 - Juin 2023	59.3	58.5	-0.1
FAPAS Food Chemistry Proficiency Test 27354 Histamine in Canned Fish - March – May 2023	120	127.5	0.8
FAPAS Food Chemistry Proficiency Test 27320 Histamine in Canned Fish - March – May 2022	168	162.1	-0.5