




Test	Assinged value (µg/kg)	Result (µg/kg)	Z-score
FAPAS Food Chemistry Proficiency Test - 17269 – Ochratoxin A in Coffee (Green) - July – August 2026	4.16	4.4	0.3
BIPEA 31B - Mycotoxins – Cereals – Rice – 7 - 4331 – June 2026	7.05	5.01	-1.01
BIPEA 31E - Mycotoxins – Feed – Soyabean Meal – 18 - 6031 – May 2026	5.07	6.70	1.07
BIPEA 31A - Mycotoxins – Dried Fruits, Spices and Other Products – Lentils – 03 - 5931 – May 2026	4.35	5.40	0.80

BIPEA 31A - Mycotoxins – Dried Fruits, Spices and Other Products – Almond Powder – 05 - 1131 – March 2026	1.38	1.50	0.29
BIPEA 31B - Mycotoxins – Cereals – Ground Corn – 19 - 0731 – December 2025	5.58	5.60	0.01
FAPAS Food Chemistry Proficiency Test - 04559 – Multi - Mycotoxins in Pet Dog Food (Dry, Cereal Based) - July – August 2025	9.33	10.7	0.7
BIPEA 31B – Mycotoxins - Cereals – Ground Rye - 4 - 5131 – June 2025	1.77	1.45	-0.6
BIPEA 31C – Mycotoxins – Baby Food - 25 - 3931 – June 2025	1.53	1.55	0.04
FAPAS Food Chemistry Proficiency Test - 04548 – Mycotoxins in Animal Feed (Cereal Based) - April – May 2025	46.8	48.0	0.1
BIPEA 31B – Mycotoxins - Cereals – Rusk - 3 - 5431 – April 2025	7.88	8.34	0.21
FAPAS Food Chemistry Proficiency Test - 04543 – Mycotoxins in Cereal Snack - February – April 2025	2.01	1.88	-0.3
BIPEA 31B – Mycotoxins - Cereals – Sorghum - 3 - 4431 – February 2025	6.00	5.77	-0.13
FAPAS Food Chemistry Proficiency Test - 04528 - Mycotoxins in Animal Feed (Cereal Based) - September – November 2024	191	170.6	-0.5

